



Does the fear of leaking exudate
mean you keep changing wound
dressings 'just in case'?

Having confidence in a dressing to manage
exudate effectively could make all the difference

www.closetozero.com

 **smith&nephew**
ALLEVYN[®] LIFE
NON-BORDERED

Foam Dressing



Supporting healthcare
professionals

ALLEVYN[◇] LIFE Non-Bordered

Foam Dressing

Unique five-layer structure is designed to:

- lock in exudate^{1,2}
- minimize the risk of skin damage^{3,4}



High fluid handling^{1,5}

- Absorbs under pressure^{1,2}
- Absorbs vertically^{*6}
- Absorbs viscous exudate⁵



Nurtures skin health

- Promotes faster healing due to the moist wound environment^{3,4}
- Minimises risk of maceration^{3,4}
- Non-adherent to the wound^{3,4}



Intelligent exudate lock in

- Locks harmful excess exudate away from the skin^{1,2}
- Helps minimize fluid leakage¹



Flexes to fit

- Cuttable and conformable^{7,8}
- Suitable for use under compression therapy⁹



Easy to apply⁹

- Soft low tack silicone aids application^{3,4}
- Minimizes wound trauma on removal^{3,4,10}

Supporting healthcare professionals for over 150 years

Smith & Nephew Pty Ltd
Australia

T +61 2 9857 3999 F +61 2 9857 3900
www.smith-nephew.com/australia

Customer Service

T 13 13 60 F 1800 671 000

Smith & Nephew Ltd
New Zealand

T +64 9 820 2840 F +64 9 820 2841
www.smith-nephew.com/new-zealand

Customer Service

T 0800 657 799 (Surgical) T 0800 807 663 (Wound)

S&N code	Size	Carton
66801748	10.5cm x 10.5cm	10 sterile dressings
66801749	16cm x 16cm	10 sterile dressings
66801751	10cm x 20cm	10 sterile dressings

References 1. Daubney L. ALLEVYN LIFE Non-Bordered - physical testing. Data on file; 2016: report DS/16/328/R v2. 2. Daubney L. Wound model of ALLEVYN Life Non-Bordered dressings. Data on file; 2016: report DS/16/470/R. 3. Rossington A, et al. Clinical performance and positive impact on patient wellbeing of ALLEVYN Life. *Wounds UK* 2013; 9: 91-95. 4. Fairlamb D. Clinical equivalence between ALLEVYN LIFE non-bordered and ALLEVYN LIFE. Data on file; 2017. 5. Daubney L. ALLEVYN LIFE Non-Bordered - fluid uptake. Data on file; 2016: report DS/16/334/Rv2. 6. Daubney L. Wound model of ALLEVYN LIFE non-bordered dressings. Data on file; 2016: report DS/16/469/R. 7. Wheldrake A. ALLEVYN Life non-bordered flexibility and conformability. Data on file; 2017: report RD/17/047. 8. Motion E. Shedding when cut - testing of Allevyn Life non-bordered and Durafiber dressings. Data on file; 2016: report DS/16/429/R. 9. Weedon F. ALLEVYN Life non-bordered claims. Data on file; 2016. 10. Davies A, et al. Selecting dressings to manage exudate and enhance patient wellbeing. *Wounds UK* 2015; 11: 54-61.

*based on performance on wound model data.

09076_AU V1 (01/18) SNI3710 ®Trademark of Smith & Nephew All Trademarks acknowledged