



## Choose the right stethoscope to help you perform at your best.

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann mechanical stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the acoustic performance, high quality and comfort you need with a Littmann stethoscope.

ACOUSTIC PERFORMANCE LEVELS						
7-8			9-10		10+	
<b>Classic II Pediatric and Infant Stethoscopes</b>			<b>Master Classic II™ Stethoscope</b>		<b>Classic III™ Monitoring Stethoscope</b>	
<b>Classic II Pediatric and Infant Stethoscopes</b>			<b>Master Classic II™ Stethoscope</b>		<b>Cardiology IV™ Diagnostic Stethoscope</b>	
<b>Classic II Pediatric and Infant Stethoscopes</b>			<b>Master Classic II™ Stethoscope</b>		<b>Master Cardiology™ Stethoscope</b>	
<b>Classic II Pediatric and Infant Stethoscopes</b>			<b>Master Classic II™ Stethoscope</b>		<b>Electronic Stethoscope Model 3200</b>	
Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.			Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.		High-performance auscultation. Hear more of what you need to hear with amplification of 24x, ambient noise reduction and optional connectivity.	
<b>Approximate Length</b>	71cm	68cm	69cm	69cm	68cm	68cm
<b>Approximate Weight</b>	105g, 95g	160g	150g	177g	175g, 185g	185g
<b>Warranty</b>	3 year	3 year	5 year	7 year	7 year	2 year
<b>Diaphragm</b>	Non-Tunable (Floating)	Tunable	Tunable (2)	Tunable (2)	Tunable	Electronic
<b>Approximate Chestpiece Diameter</b>	Double-sided Paediatric 3.4cm 2.7cm Infant 2.8cm 1.7cm	Single-sided 4.3cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm	Single-sided 4.8cm	Single-sided 5.1cm
<b>Features</b>	Specially designed for the smallest patients	Designed for easier grip and maneuvering	Convertible open bell	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Eliminates 85% (on average) of ambient noise Bluetooth® wireless technology Recording and playback



**The power to learn, in the palm of your hand**  
3M™ Littmann® Learning Institute App – Free with purchase of a 3M™ Littmann® Stethoscope!

Need more help?  
Visit [Littmann.com.au/FindMyStethoscope](http://Littmann.com.au/FindMyStethoscope)  
[Littmann.co.nz/FindMyStethoscope](http://Littmann.co.nz/FindMyStethoscope)

1300 363 878 | [www.littmann.com.au](http://www.littmann.com.au)  
0800 80 81 82 | [www.littmann.co.nz](http://www.littmann.co.nz)



# 3M™ Littmann® Stethoscope Ordering Information and Tube Colour Guide

	Classic II Infant Stethoscope	Classic II Pediatric Stethoscope	Master Classic II™ Stethoscope	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	Electronic Model 3200 Stethoscope		
All Standard Finish	2114 Black	2113 Black	2144L Black	5620 Black	5627 Burgundy	6152 Black	6156 Plum	2160 Black	3200BK27 Black
	2114R Red	2113R Red	2146 Burgundy	5621 Gray	5633 Pearl Pink	6154 Navy	6158 Raspberry	2163 Burgundy	3200BU Burgundy
	2124 Caribbean	2119 Caribbean	2147 Navy	5622 Navy	5831 Plum	6155 Hunter Green		2164 Navy	3200NB Navy
		2122 Raspberry		5623 Caribbean	5832 Lavender			2167 Plum	
				5630 Ceil Blue	5835 Turquoise				
				5839 Lemon-Lime					
All Black			2141 Black	5803 Black	5868 Burgundy	6163 Black	6168 Navy	2161 Black	
All Brass Finish									
All Copper Finish				5809 Chocolate					
All Smoke				5811 Black	5871 Raspberry	6162 Black		2176 Black	
				5812 Dark Olive		6166 Plum	6171 Turquoise		
Rainbow/Brass/Brass	2157 Raspberry	2153 Caribbean		5806 Raspberry	5807 Caribbean				
Brass/Smoke/Smoke						6164 Black			
Rainbow/Black/Black				5870 Black		6165 Black			
Mirror/Smoke/Smoke				5862 Raspberry	5863 Navy				
				5959 Ceil Blue					
Mirror/Standard/Mirror						6177 Black	6170 Burgundy		
Champagne/Smoke/Smoke				5861 Black	5864 Burgundy	6179 Black			
Champagne/Black/Champagne						6176 Burgundy	6190 Caribbean		
High-Polish Rainbow/Black/Black						6242 Navy			
High-Polish Rainbow/Smoke/Smoke						6240 Black	6241 Raspberry		
<b>COLOURED STEM SERIES: (square inset represents stem colour)</b>									
Black/Black/Colour						6200 Black	6203 Black		
						6201 Black			
High-Polish Smoke/Black/Colour						6202 Navy	6204 Black		
						6232 Black	6234 Caribbean		
High-Polish Smoke/Smoke/Colour						6238 Gray			
Rainbow/Black/Colour						6205 Plum			
Champagne/Champagne/Colour						6206 Hunter			
Smoke/Smoke/Colour				5872 Turquoise	5875 Lime				
				5873 Gray					
Mirror/Smoke/Colour				5960 Plum	5962 Pearl Pink				
Mirror/Standard/Colour				5874 Caribbean					

**Components:**  
Chestpiece/Eartube/Stem

**Standard Finish**

**Black**

**Brass Finish**

**Copper Finish**

**Smoke**

**Rainbow**

**Mirror**

**Champagne**

**High-Polish Smoke**

**High-Polish Rainbow**