

Bringing a breath of fresh air to tubular bandages

Helping you get CLOSER TO ZERO° delay in wound healing
www.closetozero.com

 **smith&nephew**

CoolGRIP™

Tubular Elastic Support Bandage

CoolFAST™

Tubular Retention Bandage



Supporting healthcare professionals



CoolGRIP™

Tubular Elastic Support Bandage

- Innovative mesh structure enhances air circulation to the skin, ensuring a more comfortable bandage, especially in warmer climates
- CoolGRIP is made from cotton and polyester, with latex free elastic threads inlaid into the bandage, providing pressure to the limb
- The fine, elastic threads are spaced close together, distributing pressure more evenly, unlike conventional bandages that incorporate thick elastic threads spaced further apart
- **Provides 14-18 mmHg sub-bandage pressure under a double layer**



Indications

- Management and treatment of lymphoedema and varicose vein conditions
- Prevention and management of gravitational oedema following cast removal
- Treatment of oedema associated with muscle sprains, strains and joint effusions
- Support of weak joints and muscles
- Retention of dressings



CoolGRIP® Tubular Elastic Support Bandage

Colour code	Product code	Product description	Limb circumference guide	Qty
A	36361527	Size A 4.5cm x 10m	Infant feet and arms	Box/1 roll
B	36361528	Size B 6.25cm x 10m	Small hands and arms	Box/1 roll
C	36361529	Size C 6.75cm x 10m	Medium arms, small ankles	Box/1 roll
D	36361530	Size D 7.75cm x 10m	Large arms, medium ankles, small knees	Box/1 roll
E	36361531	Size E 8.75cm x 10m	Large ankles, medium knees, small thighs	Box/1 roll
F	36361532	Size F 10cm x 10m	Large knees and medium thighs	Box/1 roll
G	36361533	Size G 12.5cm x 10m	Large thighs	Box/1 roll
H	36361538	Size H 14cm x 10m	Extra-large thighs	Box/1 roll

CoolFAST™

Tubular Retention Bandage

- CoolFAST has mild elastic properties provided by fine elastic threads incorporated into the fabric
- Excellent stretch and contraction properties ensure that dressings are held securely in place
- CoolFAST maintains its elasticity over a long period due to its knit structure
- CoolFAST is made from Viscose and Polyester, with latex free elastic threads inlaid into the bandage



Indications

- Dressing retention
- Protection over creams and ointments used for dermatological purposes, e.g. eczema wet wrapping



CoolFAST° Tubular Retention Bandage

Size	Product code	Product description	Limb circumference guide	Qty
RED	36361534	Red 3.5cm x 10m	10-16cm child hands and arms	Box/1 roll
GREEN	36361535	Green 5cm x 10m	14-24cm child legs and adult arms	Box/1 roll
BLUE	36361536	Blue 7.5cm x 10m	21-36cm adult arms, lower legs, knees and small thighs	Box/1 roll
Yellow	36361537	Yellow 11cm x 10m	30-50cm adult knees and thighs	Box/1 roll

Smith & Nephew Pty Ltd
Australia
T +61 2 9857 3999 F +61 2 9857 3900
www.smith-nephew.com/australia
Customer Service
T 13 13 60 F 1800 671 000

Smith & Nephew Ltd
New Zealand
T +64 9 820 2840 F +64 9 820 2841
www.smith-nephew.com/new-zealand
Customer Service
T 0800 657 799 (Surgical) T 0800 807 663 (Wound)

®All other trademarks and registered trademarks
are the property of their respective owners.

Supporting healthcare professionals for over 150 years

®Trademark of Smith & Nephew
SNI4192 (11/18)