





Littmann®
Stethoscopes



A New Class of Clarity.

The 3M™ Littmann CORE Digital Stethoscope features the outstanding acoustics, comfort and quality of a Littmann cardiology-grade stethoscope enhanced by powerful digital technology. Features like up to 40x amplification* and active noise cancellation can help you auscultate with greater confidence. The Littmann CORE Digital Stethoscope connects with Eko software on a smart device** to visualise, record and share data.

* 40x amplification based on comparison of analog versus electronically assisted amplified listening with maximum volume at the peak frequency (125 Hz).

** Smart device not included.



Connects to Eko software

Visualise, record and share heart sound waveforms by downloading the free Eko app on a smart device.**



Up to 40x amplification

Up to 40x amplification (at peak frequency, vs. analog mode).



Active noise cancellation

Noise cancellation reduces unwanted background sounds.



Analog and digital modes

Toggle between analog and amplified listening modes.



Designed for comfort

Soft-sealing ear tips provide an excellent acoustic seal and comfortable fit.



Tunable diaphragm

Tunable diaphragms on both the adult and paediatric sides of the chestpiece. The paediatric side is useful for small or thin patients, around bandages, and for carotid assessment.

Capture. Analyse. Share.

Capture, save and annotate 15, 30, 60 or 120-second recordings in the secure dashboard. The Littmann CORE Digital Stethoscope connects to Eko software on a smartphone, to create a library of heart sounds, so clinicians can monitor the progression of a disease or for education, and securely share recordings with colleagues for collaboration.



Listen first.

3M Australia Pty Limited
Building A, 1 Rivett Road
North Ryde, NSW 2113
Phone: 1300 363 878
www.littmann.com.au

3M New Zealand Limited
94 Apollo Drive, Rosedale
Auckland 0623
Phone: 0800 80 81 82
www.littmann.co.nz

Always read the label and follow the instructions.

© 3M 2020. All rights reserved.
3M, Littmann, the L Littmann logo, the L, and Cardiology IV are marks and/or registered marks of 3M Company.
Eko trademarks, including CORE, are marks and/or registered marks of Eko Devices Inc. All other trademarks are property of their respective owners.