

Budget Fridge Range

AQ MEDICAL VACCINE FRIDGES ARE A ROBUST, COST EFFECTIVE SOLUTION FOR ALL OF YOUR VACCINE REFRIGERATION NEEDS.

Key Features:

- QCPP Compliant
- Self-closing glass door with locking tab
- 24-month parts and labour warranty



AQ Medical Vaccine Fridge 145L
VMS Code: VI145MS



AQ Medical Vaccine Fridge 280L
VMS Code: VI280MS



AQ Medical Vaccine Fridge 350L
VMS Code: VI350MS

Specifications:

Model	VI145MS	VI280MS	VI350MS
Dimensions External (W x D x H) mm	595 x 595 x 820	580 x 584 x 1625	580 x 584 x 1900
Gross Litres	145L	280L	350L
Compliance	QCPP and Cold Chain compliant	QCPP and Cold Chain compliant	QCPP and Cold Chain compliant
Warranty	24 months parts and labour	24 months parts and labour	24 months parts and labour
Weight	41.5Kg	79Kg	89Kg
Operating range	2.C to 8.C	2.C to 8.C	2.C to 8.C
Shelves	4	6	6
Lock	Self closing lockable doors	Self closing lockable doors	Self closing lockable doors
Alarm	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm
Power	10 AMP	10 AMP	10 AMP
Monitoring	USB Data logger included	USB Data logger included	USB Data logger included
Cooling	Fan Assisted	Fan Assisted	Fan Assisted
Lock Included	Yes	Yes	Yes
Power Consumption	1.26KWH (Over 24hrs)	2.6KWH (Over 24hrs)	3.3KWH (Over 24hrs)