

AsGUARD[®]

PROTECTIVE

BREATHABLE

CONFORMABILITY



ACUTE
DRESSING RANGE





AsGUARD[®]

For over 30 years, Sentry Medical has been helping health care professionals provide exceptional care to their patients through a modern and quality driven portfolio.

A full range of dressings and tapes, ideal for low exuding wounds and difficult to dress areas.

AsGuard Dressings have been specifically designed to provide a secure barrier to protect against potential infection and contagions.

The extensive range comprises of various unique shapes, aiding to dress difficult areas (such as joints) and are ideal for minor wounds and first aid.

AsGUARD[®] Clear

A transparent adhesive polyurethane film with high moisture vapour permeability which adheres to the skin and forms a flexible and comfortable layer. AsGuard Clear conforms to the skin and allows easy visibility to wound healing.



AsGUARD[®] Clear+

AsGuard Clear + provides a low adherent and breathable absorbent pad permeable to moisture. The waterproof polyurethane film assists with unnecessary dressing changes and protects from infection.



AsGUARD[®] FLEX

AsGuard Flex Non-Woven Roll is a "cut to length" self-adhesive, non-woven fabric fixation tape. Highly comfortable, this tape is an ideal fixation of primary dressings for difficult contour areas without movement restriction.



AsGUARD[®] FLEX+

AsGuard Flex + Non-Woven Dressings and Rolls combine an absorbent pad with a soft and comfortable non-woven adhesive backing. The self-adhesive rolls are "cut-to-length" and ideal for fixation of difficult contour areas.



SurgiStrip[®]

SurgiStrip skin closures allow earlier post-operative suture or staple removal to assist natural healing and minimise scarring. SurgiStrips are a painless alternative for the closure of minor incisions or skin lacerations, enabling fast and easy application, with no need for anaesthetic and eliminating needle puncture marks.





AsGuard Dressing

Features & Benefits:

AsGUARD®		CLEAR DRESSING	CLEAR ROLL	CLEAR+ DRESSING	FLEX ROLL	FLEX+ DRESSING	FLEX+ ROLL	SURGI-STRIP
WATERPROOF		•	•	•				
TRANSPARENT		•	•	•				•
BREATHABLE		•	•	•	•	•	•	•
CAN BE CUT TO SIZE			•		•		•	
ABSORBENT PAD				•		•	•	
NON-WOVEN					•	•	•	•
DATE LABEL		•		•				
CONFORMABLE		•	•	•	•	•	•	•

AsGuard Dressing User Guide

AsGUARD®	CLEAR	CLEAR	CLEAR +	FLEX	FLEX +	FLEX +	SURGISTRIP
	Clear Film Dressing	Clear Film Roll	Clear Film Island Dressing	Non-Woven Roll	Non-Woven Island Dressing	Non-Woven Roll with Pad	Reinforced Skin Closures
PARTIAL THICKNESS WOUND	✓		✓		✓		
FULL THICKNESS WOUND							
LITTLE TO NO EXUDATE	✓		✓		✓		✓
MINIMAL TO LOW EXUDATE			✓		✓		✓
HIGH EXUDATE							
STAGE 1 PRESSURE ULCER	✓		✓				
STAGE 2 PRESSURE ULCER	✓		✓				
STAGE 3 PRESSURE ULCER							
REINFORCE IV SITES & LINES	✓	✓		✓		✓	
DONOR SITES	✓		✓		✓		
MILD LACERATIONS	✓		✓		✓		✓
ABRASIONS	✓		✓		✓		
SECOND DEGREE BURNS			✓		✓		
AFFIX / COVER PRIMARY DRESSING	✓	✓		✓		✓	✓
WATERPROOF BARRIER	✓	✓	✓				
NECROTIC WOUNDS	✓		✓				
PRIMARY DRESSING	✓		✓		✓		✓
SECONDARY DRESSING	✓	✓		✓		✓	
SECONDARY PADDING LAYER				✓		✓	

REORDER	DESCRIPTION	SIZE	PCS/UNIT	UNITS/CTN
ASGUARD CLEAR / FILM DRESSING				
ISB001	AsGuard Clear / Film Dressing	4cm x 5cm	50	10
ISB002	AsGuard Clear / Film Dressing	6cm x 7cm	50	10
ISB003	AsGuard Clear / Film Dressing	10cm x 12cm	25	10
ISB004	AsGuard Clear / Film Dressing	15cm x 20cm	25	10
ASGUARD CLEAR / FILM ROLL				
ISB010	AsGuard Clear / Film Roll	5cm x 10m	1	24
ISB011	AsGuard Clear / Film Roll	10cm x 10m	1	12
ASGUARD CLEAR+ / FILM ISLAND DRESSING				
ISC001	AsGuard Clear+ / Film Island Dressing	4cm x 5cm	50	10
ISC002	AsGuard Clear+ / Film Island Dressing	6cm x 7cm	50	10
ISC0915	AsGuard Clear+ / Film Island Dressing	9cm x 15cm	25	10
ISC003	AsGuard Clear+ / Film Island Dressing	10cm x 12cm	25	10
ISC028	AsGuard Clear+ / Film Island Dressing	10cm x 20cm	25	10
ISC025	AsGuard Clear+ / Film Island Dressing	10cm x 25cm	25	10
ISC026	AsGuard Clear+ / Film Island Dressing	10cm x 30cm	25	10
ISC004	AsGuard Clear+ / Film Island Dressing	15cm x 20cm	25	10
ASGUARD FLEX / NON-WOVEN ROLL				
ISD020	AsGuard Flex / Non-Woven Roll	2.5cm x 10m	1	30
ISD021	AsGuard Flex / Non-Woven Roll	5cm x 10m	1	24
ISD022	AsGuard Flex / Non-Woven Roll	7.5cm x 10m	1	24
ISD023	AsGuard Flex / Non-Woven Roll	10cm x 10m	1	12
ISD024	AsGuard Flex / Non-Woven Roll	15cm x 10m	1	12
ISD025	AsGuard Flex / Non-Woven Roll	20cm x 10m	1	6
ASGUARD FLEX+ / NON-WOVEN ISLAND DRESSING				
ISD001	AsGuard Flex+ / Non-Woven Island Dressing	7.5cm x 5cm	50	10
ISD002	AsGuard Flex+ / Non-Woven Island Dressing	8cm x 10cm	50	10
ISD030	AsGuard Flex+ / Non-Woven Island Dressing	8.25cm x 6cm	50	10
ISD003	AsGuard Flex+ / Non-Woven Island Dressing	10cm x 15cm	20	10
ISD004	AsGuard Flex+ / Non-Woven Island Dressing	10cm x 20cm	20	10
ISD006	AsGuard Flex+ / Non-Woven Island Dressing	10cm x 25cm	25	10
ISD008	AsGuard Flex+ / Non-Woven Island Dressing	10cm x 30cm	25	10
ASGUARD FLEX+ / NON-WOVEN ROLL WITH PAD				
ISD011	AsGuard Flex+ / Non-Woven Roll with Pad	5cm x 10m	1	24
ISD012	AsGuard Flex+ / Non-Woven Roll with Pad	7.5cm x 10m	1	24
ISD013	AsGuard Flex+ / Non-Woven Roll with Pad	10cm x 10m	1	12
ASGUARD FLEX+ / HEAMOPLAST-O & HEAMOPLAST				
ISD101	AsGuard Flex+ / Heamoplast-o	3.6cm Ø	100	20
ISD100	AsGuard Flex+ / Heamoplast	4cm x 5.5cm	100	20
ASGUARD SURGISTRIP / REINFORCED SKIN CLOSURES				
C-01	AsGuard SurgiStrip / Reinforced Skin Closures	3mm x 75mm	50	20
C-02	AsGuard SurgiStrip / Reinforced Skin Closures	6mm x 75mm	50	20
C-03	AsGuard SurgiStrip / Reinforced Skin Closures	6mm x 38mm	50	20
C-04	AsGuard SurgiStrip / Reinforced Skin Closures	6mm x 100mm	50	20
C-05	AsGuard SurgiStrip / Reinforced Skin Closures	13mm x 100mm	50	20
C-06	AsGuard SurgiStrip / Reinforced Skin Closures	25mm x 125mm	25	20