

Leukomed®

LEUKOMED®
WOUND DRESSINGS

SOFT, SECURE, STERILE PROTECTION



THERAPIES. HAND IN HAND.

BSN medical

LEUKOMED®

GENTLE AND ABSORBENT

Leukomed® uses Fixomull® technology, securing a highly absorbent wound pad to protect surgical and traumatic wounds.

RELIABLE AND GENTLY ADHERENT

- Provides gentle yet strong hold, protecting the wound and skin ensuring high wearing comfort so the dressing stays in place for longer
- Ensures atraumatic dressing changes assisting with patient compliance
- Rounded edges reducing roll back or snagging, allowing for undisturbed healing

HIGHLY BREATHABLE AND CONFORMING

- Soft cushioning effect for patient comfort
- Provides ideal wound healing conditions and high patient comfort assisting in longer wear time and patient comfort



SECURE ABSORPTION

- Integrated absorbent pad with low adherent contact layer absorbs exudate
- The highly absorbent wound pad allows for a reduced number of dressing changes, saving on both cost and time

STERILE AND PROTECTIVE

- Individually sealed for convenience
- Protects against friction and contaminants

Indications

- Short term use on all skin types
- Covers surgical or traumatic wounds
- Protects cuts, abrasions and lacerations

ORDERING INFORMATION

Ref-No.	Size	Pad Size	Items per Box
Leukomed®			
72380-00	5 x 7.2 cm	2.8 x 4 cm	50
72380-01	8 x 10 cm	4 x 6 cm	50
72380-02	8 x 15 cm	4 x 11 cm	50
72380-04	10 x 25 cm	5 x 20 cm	50
72380-06	10 x 35 cm	5 x 30 cm	50



BSN medical Australia & New Zealand

L6, 211 Wellington Road,
Mulgrave, Victoria, Australia 3170

www.bsnmedical.com.au
www.bsnmedical.co.nz

BSNM0004

03/2016

Customer Service Australia
T 1300 276 633 F 1300 998 830

Customer Service New Zealand
T 0508 276 111 F 0508 998 830

BSN medical