

ProShield®

Mask Range

Fluid
Resistant

PROSHIELD® CLEARVIEW

FLUID RESISTANT,
LEVEL 3 VISOR MASK

- Optimum fluid resistant
- Submicron filtration
- Protective visor
- Complies with AS4381-2015
- Suitable for all surgical procedures

THERAPIES. HAND IN HAND.

BSN medical

an Essity company

PROSHIELD® CLEARVIEW

- Wrap around visor with peripheral protection
- Submicron filtration
- Fluid Resistance (160mmHg)
- Triple pleated design for easy breathing
- Comfortable low allergy inner liner for sensitive skin
- Anti-glare strip to minimise reflection
- Suitable for use in all surgical procedures

LEVEL 3* PROTECTION

Level 3 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the maximum velocity of 160mmHg, bacterial filtration efficiency and differential pressure.

BFE

With a **Bacterial Filtration Efficiency (BFE)** of greater than 99% for media greater than 3 microns this mask provides superior protection from potential infection.

ΔP

The mask design provides excellent breathability and a **low delta P** (<2.8) ensuring coolness and comfort even when worn for an extended period.

PFE

The super high **Particulate Filtration Efficiency (PFE)** material filters >99% of particles greater than 0.1 microns. When securely fitted, the mask helps to provide an effective barrier to particulate matter.



The **Fluid Resistant** qualities of the mask provide protection against fluid (160mmHg/21.3kPa) strikethrough.

APPLICATION

For all surgical procedures, major trauma first aid or in any area where the health care worker is at risk of blood or bodily fluid splash (e.g. orthopaedic, cardiovascular procedures).



FLUID RESISTANCE

- The fluid resistance of a surgical face mask is the ability of the mask to limit the passage of blood or body fluids from the outer facing through the mask layer to the inner facing
- The pressure of the spurt is dependent on a number of factors. These include:
 - ▶ Blood pressure of the patient
 - ▶ Velocity of the spurt
 - ▶ Distance between the source of the spurt to the mask wearer

ORDERING INFORMATION

Ref-no.	Description	Items per unit	Shipper Qty
72509-11	Clearview Visor Masks	bx/50	6 bxs

*AS 4381:2015 Single-use face masks for use in healthcare

BSN medical Australia & New Zealand

L6, 211 Wellington Road,
Mulgrave, Victoria, Australia 3170

www.bsnmedical.com.au
www.bsnmedical.co.nz

BSN0343

09/2017

Customer Service Australia
T 1300 276 633 F 1300 998 830

Customer Service New Zealand
T 0508 276 111 F 0508 998 830

BSNmedical

an Essity company