

# ProShield®

Mask Range

Fluid  
Resistant

## PROSHIELD® DUCKBILL FR

SINGLE PLEATED, LEVEL 2  
SUBMICRON FILTRATION  
FLUID RESISTANT MASK

- Submicron filtration
- Fluid resistant
- Excellent breathability
- Complies with AS4381: 2015

THERAPIES. HAND IN HAND.

**BSN** medical

an Essity company

# PROSHIELD® DUCKBILL FLUID RESISTANT MASK

- Submicron filtration
- Fluid Resistance (120mmHg)
- Comfortable low allergy inner liner for sensitive skin
- Single pleated Duckbill design for easy breathing
- Soft and durable ties

## LEVEL 2\* PROTECTION

Level 2 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the medium velocity of 120mmHg, bacterial filtration efficiency and differential pressure.

**BFE**

With a **Bacterial Filtration Efficiency (BFE)** of 98% for media greater than 3 microns, this mask provides superior protection from potential infection.

**▲ P**

The mask design provides excellent breathability and a **low delta P** ( $\leq 1.5$ ) ensuring coolness and comfort even when worn for an extended period.

**PFE**

The super high **Particulate Filtration Efficiency (PFE)** material filters  $\geq 98\%$  of particles greater than 0.1 microns. When securely fitted, the mask helps to provide an effective barrier to particulate matter.



The **Fluid Resistant** qualities of the mask provide protection against fluid (120mmHg/16.0kPa) strikethrough.

## APPLICATION

For use in emergency departments, dentistry, changing dressings on small wounds or healing wounds and where minimal blood droplet exposure may possible occur (e.g. endoscopy procedures).



## FLUID RESISTANCE

- The fluid resistance of a surgical face mask is the ability of the mask to limit the passage of blood or body fluids from the outer facing through the mask layer to the inner facing.
- The pressure of the spurt is dependent on a number of factors. These include:
  - ▶ Blood pressure of the patient
  - ▶ Velocity of the spurt
  - ▶ Distance between the source of the spurt to the mask wearer

## ORDERING INFORMATION

Ref-no.	Description	Items per unit	Shipper Qty
72509-00	ProShield® Duckbill FR	bx/50	30 bxs

\*AS 4381:2015 Single-use face masks for use in healthcare

### BSN medical Australia & New Zealand

L6, 211 Wellington Road,  
Mulgrave, Victoria, Australia 3170

[www.bsnmedical.com.au](http://www.bsnmedical.com.au)  
[www.bsnmedical.co.nz](http://www.bsnmedical.co.nz)

BSN0348

10/2017

Customer Service Australia  
T 1300 276 633 F 1300 998 830

Customer Service New Zealand  
T 0508 276 111 F 0508 998 830

**BSN** medical

an Essity company