



# Optimal Refrigeration Range

# Premium Fridge Range

REMOVE THE RISK OF STAFF INCORRECTLY DOWNLOADING DATA VIA A MANUAL DATA LOGGER WITH AUTOMATIC LOGGING VIA THE MEDISAFE PLUS BUILT IN CLOUD MONITORING SYSTEM.

## Key Features:

- 12 Months free cloud based monitoring
- View and Print temperature data
- Text and Email alerts
- Built in AQ box



Medisafe Plus 80  
VMS Code: VI80MS



Medisafe Plus 311  
VMS Code: VI311MS



Medisafe Plus 371  
VMS Code: VI371MS

## Specifications:

Model	VI80MS	VI311MS	VI371MS
<b>Dimensions External (W x D x H) mm</b>	500 x 560 x 830	595 x 595 x 1630	595 x 595 x 1930
<b>Gross Litres</b>	80L	306L	381L
<b>Compliance</b>	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested
<b>Warranty</b>	24 months parts and labour	24 months parts and labour	24 months parts and labour
<b>Weight</b>	42Kg	78Kg	89Kg
<b>Operating range</b>	2.C to 8.C	2.C to 8.C	2.C to 8.C
<b>Shelves</b>	3	5	6
<b>Lock</b>	Self closing lockable doors	Self closing lockable doors	Self closing lockable doors
<b>Alarm</b>	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm
<b>Power</b>	10 AMP	10 AMP	10 AMP
<b>Monitoring</b>	AQ Box Cloud based monitoring (initial 12 months free)	AQ Box Cloud based monitoring (initial 12 months free)	AQ Box Cloud based monitoring (initial 12 months free)
<b>Cooling</b>	Fan Assisted + Door Fan Switch	Fan Assisted	Fan Assisted
<b>Body Type</b>	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External
<b>Power Consumption</b>	3.08 KWH (over 24hrs)	3.69KWH (over 24hrs)	4.71 KWH (over 24hrs)

# Premium Fridge Range

**EXCEPTIONAL QUALITY  
AND TRULY OUTSTANDING  
FEATURES, AT AN  
AFFORDABLE PRICE.**

## Key Features:

- QCPP Compliant
- Independently tested by The Pharmacy Guild of NSW
- Touch Screen Digital Controller
- Audible & Visual Alarm
- Built In Lock
- High Ambient Rating 40°C / 75% Humidity \*
- 3-year comprehensive warranty \*
- VS1300 Premium different specs



Vacc-Safe 250LT Fridge  
VMS Code: VACCSAFE2



Vacc-Safe 400LT Fridge  
VMS Code: VACCSAFE4



Vacc-Safe 600LT Fridge  
VMS Code: VACCSAFE6



Vacc-Safe 1000LT Fridge  
VMS Code: VACCSAFE10



Vacc-Safe 1300LT Fridge  
VMS Code: VACCSAFE13

## Specifications:

Model	VACCSAFE2	VACCSAFE4	VACCSAFE6	VACCSAFE10	VACCSAFE13
<b>Dimensions External (W x D x H) mm</b>	560 x 640 x 1454	560 x 640 x 1930	700 x 710 x 1930	1050 x 710 x 1950	1314 x 845 x 2105
<b>Gross Litres</b>	225L	400L	Gross: 600L Net: 500L	1000L	1300L
<b>Compliance</b>	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested
<b>Warranty</b>	36 months parts and labour + 60 months parts only on compressor	36 months parts and labour + 60 months parts only on compressor	36 months parts and labour + 60 months parts only on compressor	36 months parts and labour + 60 months parts only on compressor	24 months parts and labour
<b>Weight</b>	65Kg	87Kg	110Kg	130Kg	190Kg
<b>Operating range</b>	2.C to 8.C	2.C to 8.C	2.C to 8.C	2.C to 8.C	2.C to 8.C
<b>Shelves</b>	5	7	6	12 (6 per door)	12 (6 per door)
<b>Lock</b>	Self closing lockable doors	Self closing lockable doors	Self closing lockable doors	Self closing lockable doors	Self closing lockable doors
<b>External Alarm Contacts</b>	NO/NC External Alarm Contacts Fitted	NO/NC External Alarm Contacts Fitted	NO/NC External Alarm Contacts Fitted	NO/NC External Alarm Contacts Fitted	NO/NC External Alarm Contacts Fitted
<b>Alarm</b>	High / Low & Door Open Audible & Visual Alarm	High / Low & Door Open Audible & Visual Alarm	High / Low & Door Open Audible & Visual Alarm	High / Low & Door Open Audible & Visual Alarm	High / Low & Door Open Audible & Visual Alarm
<b>Power</b>	10 AMP	10 AMP	10 AMP	10 AMP	10 AMP
<b>Monitoring</b>	USB Log Tag data logger + free software download	USB Log Tag data logger + free software download	USB Log Tag data logger + free software download	USB Log Tag data logger + free software download	USB Log Tag data logger + free software download
<b>Cooling</b>	Fan Assisted	Fan Assisted	Fan Assisted	Fan Assisted	Ventilated
<b>Body Type</b>	White Powder Coated Steel Internal & External	White Powder Coated Steel Internal & External	White Powder Coated Steel Internal & External	White Powder Coated Steel Internal & External	Stainless Steel Exterior / Interior
<b>Power Consumption</b>	1.428KWH (over 24hrs)	3.92KWH (over 24hrs)	3.7KWH (over 24hrs)	7.4KWH (over 24hrs)	2.3KWH (Over 24hrs)

# Intermediate Fridge Range

COMBINING SUPERB  
PERFORMANCE AND QUALITY  
WITH GREAT VALUE.

## Key Features:

- QCPP Compliant
- Independently tested by the Pharmacy Guild Of NSW
- Self-closing glass door with inbuilt lock
- 24-month comprehensive warranty
- Rated to 38 Degrees



Medi Guard® 401 Fridge  
VMS Code: VACCMED4



Medi Guard® 601 Fridge  
VMS Code: VACCMED6



Medi Guard® 1001 Fridge  
VMS Code: VACCMED10

## Specifications:

Model	VACCMED4	VACCMED6	VACCMED10
<b>Dimensions External (W x D x H) mm</b>	606 x 620 x 1940	700 x 790 x 2058	11130 x 720 x 2023
<b>Gross Litres</b>	380L	600L	1000L
<b>Compliance</b>	QCPP and Cold Chain compliant + Pharmacy Guild Tested	TBC	QCPP and Cold Chain compliant + Pharmacy Guild Tested
<b>Warranty</b>	24 months parts and labour	24 months parts and labour	24 months parts and labour
<b>Weight</b>	75Kg	110Kg	142Kg
<b>Operating range</b>	2.C to 8.C	2.C to 8.C	2.C to 8.C
<b>Shelves</b>	7	6	12
<b>Lock</b>	Self closing lockable doors	Self closing lockable doors	Self closing lockable doors
<b>Alarm</b>	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm
<b>Power</b>	10 AMP	10 AMP	10 AMP
<b>Monitoring</b>	Medi Log 2 Data Logger Audible and Visible Alerts	Medi Log 2 Data Logger Audible and Visible Alerts	Medi Log 2 Data Logger Audible and Visible Alerts
<b>Cooling</b>	Fan Assisted	Fan Assisted	Fan Assisted
<b>Body Type</b>	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External
<b>Power Consumption</b>	4.935KWH (Over 24hrs)	6.209KWH (Over 24hrs)	9.903KWH (Over 24hrs)

# Budget Fridge Range

**AQ MEDICAL VACCINE FRIDGES ARE A ROBUST, COST EFFECTIVE SOLUTION FOR ALL OF YOUR VACCINE REFRIGERATION NEEDS.**

## Key Features:

- QCPP Compliant
- Self-closing glass door with locking tab
- 24-month parts and labour warranty



AQ Medical Vaccine Fridge 145L  
VMS Code: VI145MS



AQ Medical Vaccine Fridge 280L  
VMS Code: VI280MS



AQ Medical Vaccine Fridge 350L  
VMS Code: VI350MS

## Specifications:

Model	VI145MS	VI280MS	VI350MS
<b>Dimensions External (W x D x H) mm</b>	595 x 595 x 820	580 x 584 x 1625	580 x 584 x 1900
<b>Gross Litres</b>	145L	280L	350L
<b>Compliance</b>	QCPP and Cold Chain compliant	QCPP and Cold Chain compliant	QCPP and Cold Chain compliant
<b>Warranty</b>	24 months parts and labour	24 months parts and labour	24 months parts and labour
<b>Weight</b>	41.5Kg	79Kg	89Kg
<b>Operating range</b>	2.C to 8.C	2.C to 8.C	2.C to 8.C
<b>Shelves</b>	4	6	6
<b>Lock</b>	Self closing lockable doors	Self closing lockable doors	Self closing lockable doors
<b>Alarm</b>	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm
<b>Power</b>	10 AMP	10 AMP	10 AMP
<b>Monitoring</b>	USB Data logger included	USB Data logger included	USB Data logger included
<b>Cooling</b>	Fan Assisted	Fan Assisted	Fan Assisted
<b>Lock Included</b>	Yes	Yes	Yes
<b>Power Consumption</b>	1.26KWH (Over 24hrs)	2.6KWH (Over 24hrs)	3.3KWH (Over 24hrs)







**Mail:** PO Box 100, Kingsgrove, NSW 1480

**Web:** [www.vitalmed.com.au](http://www.vitalmed.com.au)

**Email:** [orders@vitalmed.com.au](mailto:orders@vitalmed.com.au)

**Phone:** 1300 557 651

**Fax:** 1300 557 631

**Online Ordering:** [www.vitalmedicalsupplies.com.au](http://www.vitalmedicalsupplies.com.au)